


PiC Diabetes
care



Insupen

Every day is problem-free.



Pic Solution.
Specialists in self care.

*PiC Diabetes Care:
as always, your key to freedom.*

PiC has always been close to people. Its 60 years experience in selfcare is based on a thorough knowledge of customers and on the ability to understand their deepest needs.



This attention and ongoing **research** have created **PiC Diabetes Care**, a complete range of simple and innovative products **designed to help manage diabetes** and its complications **freely**, at any time.

PiC Diabetes Care means **scientific research, sensitivity and experience.**



MADE IN ITALY





Insupen pen needles: maximum performance, minimum pain.

PiC **Insupen** is the pen needles range that matches the reduction of the needle outer diameter **with a broad variety** of lengths in order to meet **almost all injection needs**. **Insupen pen needles are the result of the indolor™ experience**, a legacy of research and expertise that guides PiC Diabetes Care in developing an **indolor™ therapy**, a goal achieved thanks to many breakthroughs.

- 1 Insupen32G** *a needle with a diameter of 0,23mm available in 3 lengths (4, 6 and 8 mm).*
- 2 Insupen33G** *the diameter of 0.20 mm is the result of continuous technological development. Clinically tested.**
- 3 EXTR 3 ME** *the exclusive needle geometry that reduces the penetration force needed for the injection.*
- 4 HIGH FLOW THIN WALL** *an even wider lumen for an enhanced flow (vs. Insupen thin wall needles).*

*Data on file.

1. Insupen 32G: the complete line that caters to all habits.

When talking about injections, habits make a huge difference. Thanks to its broad range (4, 6 and 8 mm), **Insupen 32G** is the 32-gauge pen needle **able to meet almost all needs**, offering the **freedom to choose an indolor™ injection without changing** the injection habit.



Thin needle available in 3 lengths

More comfort, wider choice

Pic Solution.
Specialists in self care.

2. Insupen 33G: the thinnest Insupen pen needle.




Minimal dimensions

Our research brings you **Insupen 33G**, the cutting edge of **inZolor®** technology and **comfort**: its diameter of just 0,20 mm and length of 4 mm **allow for an inZolor® injection** and comfortable therapy.

Maximum adaptability

Insupen 33G offer most patients an optimal therapy option. The average dermis thickness is between 1,8 and 2,5 mm, therefore a **4 mm needle** will reach the right administration point in most patients.



Advanced technology

The small size, the new Extr3me tip geometry and the High Flow Thin Wall technology make **the 33G pen needle a genuine key to freedom** for diabetics. Freedom from the fear of injections, for **a comfortable therapy**.





3. New Extr3me needle geometry.

Insupen presents the new range of pen needles developed from the **indolor**® experience **research and featuring Extr3me's exclusive three cutting edges**, available on 33G, 32G and 31G.

Insupen

EXTR 3 ME

Less invasiveness


The **revolutionary design** of Insupen Extreme pen needles **significantly reduces the needle size** and hence the feeling of **invasiveness and discomfort**.



Comparison Insupen Extr3me vs Insupen Traditional.
Cross-section of less material in red.

Less penetration force

Extreme's **sharper bevels** make the needle penetrate more easily into the skin and with less force **thus causing less pain**.



Less invasiveness
Less penetration force*

Insupen 33G, 32G and 31G adopt the new Extreme needle geometry.

*Comparison Insupen Extr3me vs Insupen Traditional.


Pic Solution.
Specialists in self care.

4. Enhanced flow.

Thin walled needles feature a **wider internal diameter**, that allow to **reduce the insulin administration time**.

High Flow Thin Wall: the details that make a difference

With an inner wall of just 0,045 mm, and a greater internal lumen, the High Flow Thin Wall technology **provides an enhanced drug flow** thus **reducing administration time**, without the need of increasing the needle outer diameter.



*Insupen 32Gx4 Extreme vs Insupen 32Gx4 traditional. Data on file.

**High Flow Thin Wall Technology:
wider inner gauge**

Shorter administration time**

**Insupen 32Gx4 Extreme vs Insupen 32Gx4 traditional. Data on file.



**5. Approved by those
who tried it.**

Extr3me keeps promises of increasingly comfortable injections: a **genuine indolor[®]** experience.



Pic Solution.
Specialists in self care.

Insupen pen needles: a comprehensive range for any need.

Insupen33G

EXTR 3.ME

33G
4 mm



Insupen32G

EXTR 3.ME



Insupen31G

EXTR 3.ME



Insupen30G



63 099621 080 020

108866.M1.1

www.picsolution.com

Manufactured by: Distributed by:



Artsana S.p.A.
Via Saldarini Catelli, 1
22070 Grandate (CO) Italy
Phone: +39 031-382.111
Fax: +39 031-382.400
pic.solution@salessupport.com

DOMREX PHARMA INC.
3232 Ave. Jacques-Bureau
Laval (Quebec) Canada H7P0A9
Tel: 450-687-0990,
1-888-687-0990
Email: info@domrexpharma.com
www.domrexpharma.com



Specialists in self care.